




Novel Diagnostics for improved colorectal cancer staging



**Reaching beyond histopathology
by using molecular biomarkers**

With current diagnostics, one fourth of colorectal cancer patients judged to be cured by surgery will die from tumor recurrence

How can we reduce recurrence rates and increase survival?

Lymph nodes are key elements of the TNM (Tumor, Node and Metastasis) staging System. The presence or absence of disseminated tumor cells is the single most important factor when predicting disease-free survival and overall survival for patients with colorectal cancer without distant metastasis. The lymph node status is also a crucial factor when deciding on adjuvant chemotherapy treatment after curative surgery.

However, today only about 1 % of each lymph node is analyzed in clinical routine.

ColoNode® is a novel In Vitro Diagnostics for colorectal cancer that detects and characterizes tumor cells in lymph nodes by measuring mRNA levels of five biomarkers. ColoNode® estimates the risk of tumor recurrence based on the expression profile of these biomarkers. ColoNode® can analyze up to 100% of the lymph node volume. This provides an increased basis for decision-making in staging and decisions on adjuvant chemotherapy treatment and follow-up.

Today only about

1 % of each lymph node is analyzed in clinical routine





ColoNode®

ColoNode® is a multiplex qRT-PCR assay that analyzes the mRNA levels of five biomarkers and one control gene:

- CEACAM5 protein is a well-established tumor marker in colorectal cancer. CEACAM5 mRNA levels correlate strongly with the number of tumor cells^{1,2}
- Expression of KLK6 mRNA is associated with poor prognosis³
- Expression of SLC35D3 mRNA is associated with poor prognosis⁴
- High levels of stromal POSTN mRNA correlate to poor prognosis⁴
- Lymph nodes harboring tumor cells with high levels of MUC2 mRNA is an indicator of good prognosis²
- ColoNode gives an estimate of risk of recurrence and allocates the patient to one of four risk groups⁵

Most important scientific results



ColoNode is superior to histopathology for detection of tumor cells in lymph nodes

Histopathology (H&E) and CEACAM5 mRNA level by ColoNode were compared for detection of tumor cells in 2393 lymph nodes of 196 colon cancer patients.

ColoNode **increased detection of metastases and micrometastases 3.34-fold** compared to histopathology.

73% of the lymph nodes were CEACAM5-positive, that is, harbored tumor cells.

References 5 and 6


Allocation to risk groups by ColoNode

For each CEACAM5-positive lymph node a risk factor is calculated by a formula for the expression of biomarkers KLK6, SLC35D3, POSTN and MUC2.

(-1) and (0) = low risk (1) = medium risk
(2) = high risk (3) = very high risk


The lymph node with the highest risk factor defines the risk group to which the patient is allocated

Risk factors of individual lymph nodes of a colon cancer patient



ColoNode surpasses histopathology in identifying patients at risk of recurrence

The COLONODE-study – a prospective multicenter study in which all lymph nodes examined by histopathology in the clinic also were analyzed by ColoNode. Lymph nodes were from 196 colon cancer patients receiving curative surgery. At 3-year follow-up 36 patients had died from their cancer or lived with recurrent disease.



COLONODE RISK GROUP 1 OR HIGHER (1,2,3):

- identified **all** patients who had recurred that were identified by histopathology (21 patients; 58%)
- identified **9 additional** patients who had recurred (**25%, P < 0.0001**)
- proved to be a highly significant, independent risk factor with a **hazard ratio of 4.24** (95% CI, 1.42-12.69; P = 0.01) in **multivariate analysis**, while histopathology (pTN-stage III vs I/II) lost its univariate significance

PATIENTS IN COLONODE RISK GROUP LOW (0,-1), = 51%:

- have very low risk of recurrence, that is, these patients **are likely cured by surgery alone**
- PATIENTS IN COLONODE RISK GROUP HIGH/VERY HIGH (2,3), = 29%:
- are at high risk of recurrence with a hazard ratio of 10.6 (95% CI: 3.6-31.2) compared to the Low (0,-1) group. These patients **are likely to benefit from postoperative adjuvant treatment**

Reference 6

References

1. Öberg Å, Lindmark G, Israelsson A, Hammarström S, Hammarström M-L. Detection of occult tumor cells in lymph nodes of colorectal cancer patients using real-time quantitative RT-PCR for CEA and CK20 mRNAs. *Int J Cancer*. 2004;111:101-110. doi: 10.1002/ijc.20231. PMID: 15185350.
2. Ohlsson L, Israelsson A, Öberg Å, Palmqvist R, Stenlund H, Hammarström M-L, Hammarström S, Lindmark G. Lymph node CEA and MUC2 mRNA as useful predictors of outcome in colorectal cancer. *Int J Cancer*. 2012;130:1833-1843. doi: 10.1002/ijc.26182. PMID: 21618511.
3. Ohlsson L, Lindmark G, Israelsson A, Palmqvist R, Öberg Å, Hammarström M-L, Hammarström S. Lymph node tissue kallikrein-related peptidase 6 mRNA – a progression marker for colorectal cancer. *Br J Cancer*. 2012;107:150-157. doi: 10.1038/bjc.2012.220. PMID: 22699826.
4. Olsson L, Hammarström M-L, Israelsson A, Lindmark G, Hammarström S. Allocating colorectal cancer patients to different risk categories by using a five-biomarker mRNA combination in lymph node analysis. *PLoS ONE*. 2020;15:e0229007. doi: 10.1371/journal.pone.0229007. eCollection 2020. PMID: 32049988.
5. Olsson LM, Lindmark GE, Israelsson ACE, Korkocic D, Hammarström SG, Hammarström M-LKC. CEACAM5, KLK6, SLC35D3, POSTN and MUC2 mRNA analysis improves detection and allows characterization of tumor cells in lymph nodes of colon cancer patients. *Dis Colon Rectum*. 2021;64:1354-1363. doi: 10.1097/DCR.0000000000002151. PMID: 34192710.
6. Lindmark G, Olsson L, Sitohy B, Israelsson A, Blomqvist J, Kero S, Roshdy T, Söderholm M, Turi A, Isaksson J, Sakari T, Dooper M, Dafnis G, Forsberg P, Skovsted S, Walldén M, Kung CH, Rutegård M, Nordmyr J, Muhrbeck M, Hammarström S, Hammarström ML. qRT-PCR analysis of CEACAM5, KLK6, SLC35D3, MUC2 and POSTN in colon cancer lymph nodes-An improved method for assessment of tumor stage and prognosis. *Int J Cancer*. 2024;154: 573-584. doi: 10.1002/ijc.34718. Epub 2023 Sep 12. PMID: 37700602.



1.



2.



3.




4.



5.



6.



Distributed by
BIOCARTIS

Distributed in the EU by:

Biocartis NV | Generaal De Wittelaan 11B | 2800 Mechelen
Belgium | Phone: +3215632600 | biocartis.com



Publications

www.colonode.se/publications

www.colonode.se

HILOPROBE

HiloProbe AB, Tvistevägen 48C, 907 36 Umeå, Sweden

www.hiloprobe.se