

IDYLLA™ SPECIMEN REQUIREMENTS



One of the biggest challenges in oncology biomarker testing is the ability to obtain samples of sufficient size and quality. The Idylla™ System only needs a minimal amount of sample.

SOLID BIOPSY SPECIMEN REQUIREMENTS

IVD Tests



MSI 510(k)¹

1 x 4 µm FFPE tissue section → 62.5-750 mm²
 1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 33%, if less, macrodissection is required

RUO Assays

EGFR

1 x 5 µm FFPE tissue section
 Neoplastic cells ≥ 10% - if less, macrodissection is required

GeneFusion

1 x 5 µm FFPE tissue section if tissue area ≥ 20 mm²
 3 x 5 µm FFPE tissue section if tissue area < 20 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

BRAF

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 50% - if less, macrodissection is required

KRAS

NRAS-BRAF-EGFR S492R

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

POLE-POLD1

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

IDH1-2

Extracted DNA: 50 µl (concentration ≥ 10 ng/µl)
 FFPE: maximum 3 FFPE tissue sections
 Neoplastic cells ≥ 10% - if less, macrodissection is required
 Whole Blood or Bone Marrow: 10 µl

PIK3CA-AKT1

1 x 5 µm FFPE tissue section if tissue area 50-600 mm²
 1 x 10 µm FFPE tissue section if tissue area 25-300 mm²
 Neoplastic cells ≥ 20% - if less, macrodissection is required

LIQUID BIOPSY SPECIMEN REQUIREMENTS²

RUO Assays



ctEGFR

2 ml plasma collected using EDTA tubes

ctKRAS

ctNRAS-BRAF-EGFR S492R

1 ml plasma collected using either EDTA or Streck tubes

(1) For in vitro diagnostic use. For use on the Biocartis Idylla™ System only. The Idylla™ MSI Test, for use on the Idylla™ System, uses formalin-fixed, paraffin-embedded (FFPE) tissue sections of human CRC tumor, from which nucleic acids are liberated, then analyzed using PCR amplification of seven monomorphic biomarkers (ACVR2A, BTBD7, D1D01, MRE11, RYR3, SEC31A and SULF2) and subsequent melt-curve analysis. The Idylla™ MSI Test reports results as either microsatellite stable (MSS), or microsatellite instability high (MSI-H) or invalid. The Idylla™ MSI Test is indicated for use by healthcare professionals for the qualitative identification of microsatellite instability (MSI) in colorectal cancer (CRC) tumors, indicative of mismatch repair deficiency, and as an aid in the identification of probable Lynch syndrome to help identify patients that would benefit from additional genetic testing to diagnose Lynch syndrome. The results from the Idylla™ MSI Test should be interpreted by healthcare professionals in conjunction with other clinical findings, family history, and other laboratory data. The Idylla™ MSI Test should not be used for diagnosis of CRC. The clinical performance of this device to guide treatment decision for MSI high patients has not been established.

(2) K₂EDTA tubes: process within 4 hours. Streck Cell-Free DNA BCT® Tubes: process within 3 days.



Please refer to the product instructions for complete instructions on how to process samples on Idylla™.

Idylla™ Platform is listed as a class II device in the US under establishment registration 3009972873. Idylla™ MSI Test is cleared in the US under K211181 for Lynch Syndrome. Idylla™ EGFR, ctEGFR, BRAF, KRAS, ctKRAS, NRAS-BRAF-EGFR S492R, ctNRAS-BRAF-EGFR S492R, POLE-POLD1 & PIK3CA-AKT1 Mutation Assays; Idylla™ GeneFusion Assay; and Idylla™ IDH1-2 Mutation Assay Kit are for Research Use Only, not for use in diagnostic procedures. For more information on the acceptable use and licenses of Idylla™ products, please visit www.biocartis.com/en/license-statements.

Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis NV. Idylla™ is available for sale in Europe, the US and many other countries. Please check availability with a Biocartis representative.

© March 2025, Biocartis NV. All rights reserved.