

BECAUSE  
TIME MATTERS

# IDYLLA™ POLE-POLD1 MUTATION ASSAY

For the detection of the  
hypermuted phenotype  
associated with pathogenic  
mutations in POLE and POLD1



THINK IDYLLA™  
BECAUSE TIME MATTERS

# INTRODUCING THE IDYLLA™ POLE-POLD1 MUTATION ASSAY!



Qualitative detection of **17 POLE** pathogenic mutations across 4 exons:

- Exon 9: P286H, P286L, P286R, P286S, M295R, S297F
- Exon 11: F367S, F367V, D368Y
- Exon 13: V411L (G>C; G>T), L424I, L424V, P436R, M444K
- Exon 14: A456P, S459F

Additionally, **1 POLD1** mutation is detected in exon 12: S478N



**Fully automated** molecular testing platform

**On-demand** testing, without the need for sample batching



**Under 3 minutes** hands-on time

Assay turnaround time (TAT) of approx. **95 minutes**



Directly from 1 **FFPE** tissue section


## SPECIMEN REQUIREMENTS

- 50-600mm<sup>2</sup> tissue area for 5 µm FFPE tissue sections
- 25-300mm<sup>2</sup> tissue area for 10 µm FFPE tissue sections
- ≥ 10% neoplastic cells



# A MULTI-CENTER STUDY HIGHLIGHTS THE ASSAY'S POTENTIAL IN ENDOMETRIAL CANCER RESEARCH

A retrospective, multi-center study was conducted to evaluate the Assay's potential by comparing it to NGS, Sanger of qPCR using a cohort of endometrial cancer samples<sup>1</sup> (Barault et al., 2024).



1: MSKCC (USA); 2: Compunet (USA); 3: CMOCO (USA); 4: Methodist (USA); 5: Augusta (USA); 6: Hvidovre (DK); 7: Arnau de Vilanova (ES); 8: Ludwigsburg (DE); 9: Kassel (DE); 10: Graz (AT)

## IDYLLA™ DEMONSTRATED 98.6% ACCURACY WITH ENDOMETRIAL CANCER SAMPLES

434 tissue samples met all inclusion criteria<sup>2</sup> and generated valid results with Idylla™ and the reference methods.

**Table 1.** Aggregated results comparing Idylla™ and the reference methods.

		Reference methods		
		Mutated	Wild Type	Total
Idylla™	Mutated	175	2*	177
	Wild Type	4**	253	257
	Total	179	255	434

\* 2/2 false positive results were potentially due to the use of a lower sensitivity reference method (Sanger).

\*\* 4/4 false negative results were obtained with a low amount of amplifiable DNA.


**Table 2.** Concordance rates of Idylla™ with the reference methods.

	All samples (n=434)	Without low input samples (n=386)
PPA	97.7%	100%
NPA	99.2%	99.2%
OPA	98.6%	99.5%

PPA: positive percent agreement, NPA: negative percent agreement, OPA: overall percent agreement.

(1) The assay used during this prototype study included 12 POLE mutations. Five POLE mutations (L424V, F367V and P286H/L/S) as well as the POLD1 mutation (S478N) were added after data collection to increase mutation coverage.

(2) ≥ 10% tumor cells, ≥ 0.25 mm<sup>3</sup> tissue area.



**Figure 1.** Distribution of the detected mutations for the 175 mutated samples<sup>3</sup>.

## CONCLUSION



Fully automated solution to detect 99% of the known pathogenic POLE and POLD1 mutations



Demonstrated 98.6% accuracy with 434 endometrial cancer tissue samples

The Idylla™ POLE-POLD1 Mutation Assay provides a rapid and reliable solution to close an important gap in oncology research by enabling molecular classification of endometrial cancer samples.

## REFERENCES

Barault, L. et al. (2024). The Idylla™ POLE Mutation Assay, A New Tool For Direct Mutation Detection From FFPE Tissue. *AMP 2024 Annual Meeting*.

Biocartis NV  
 Generaal De Wittelaan 11B  
 2800 Mechelen - Belgium  
 +32 15 632 888

Biocartis US Inc.  
 2 Pierce Place, Suite 1510  
 Itasca, IL 60143 - US  
 +1 (844) 443-9552

Follow us on     
[www.biocartis.com](http://www.biocartis.com)  
[customerservice@biocartis.com](mailto:customerservice@biocartis.com)  
[customerserviceUS@biocartis.com](mailto:customerserviceUS@biocartis.com)

(3) Idylla™ reports the different mutations per exon.

Idylla™ POLE-POLD1 Mutation Assay is for Research Use Only (RUO), not for use in diagnostic procedures. Idylla™ Platform is CE-marked in Europe in compliance with EU IVD Regulation 2017/746, listed as a class II device in the US under establishment registration 3009972873, and registered in many other countries. Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis NV. Idylla™ is available for sale in Europe, the US and many other countries. Please check availability with a Biocartis representative. ©March 2025, Biocartis NV. All rights reserved.